



IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

Attorney Docket No.: P-576 (TI-0022)
Inventors: Huber et al.
Serial No.: 09/770,410
Filing Date: January 25, 2001
Examiner: Therkorn, Ernest G.
Customer No.: 26259
Group Art. Unit: 1723
Confirmation No.: 6186
Title: Method and Apparatus for Separating
Polynucleotides Using Monolithic
Capillary Columns

"Express Mail" Label No. EV735453368US

Date of Deposit: March 6, 2006

I hereby certify that this paper is being deposited
with the United States Postal Service "Express Mail
Post Office to Addressee" service under 37 CFR 1.10
on the date indicated above and is addressed to the
Commissioner for Patents, Mail Stop AF,
P. O. Box 1450, Alexandria, VA 22313-1450.

By Jane Massey Licati
Typed Name: Jane Massey Licati, Reg. No. 32,257

Commissioner for Patents
Mail Stop AF
P.O. Box 1450
Alexandria, VA 22313-1450

Dear Sir:

DECLARATION UNDER RULE 9 1.131

I, Andreas Premstaller, hereby declare that:

1. I am a co-inventor, together with Christian Huber
and Herbert Oberacher, in U.S. Patent Application Serial No.
09/770,410 filed June 7, 2000 and am most familiar with the
subject matter of this application and the research effort which
lead to the discovery of the instant invention. All the work
described in the following paragraph occurred at the Institute of

BEST AVAILABLE COPY

Attorney Docket No.: P-576 (TI-0022)
Inventors: Huber et al.
Serial No.: 09/770,410
Filing Date: January 25, 2001
Page 2

Analytical Chemistry and Radiochemistry in Innsbruck, Austria, a recognized WTO member country since January 1, 1995.

2. I have reviewed Gusev et al. ((September 1999) *J. Chromatography* 855:273-293) and find that this reference describes a porous monolithic packing prepared with polystyrene-divinylbenzene support which is covalently attached to a fused silica capillary inner wall treated with a coupling agent trimethoxysilyl propyl methacrylate to provide anchoring sites for grafting of the polymer to the silica surface. The median pore radius for a monolithic sample prepared with ethanol is, as estimated by Gusev, about 5 micrometers.

3. Our invention referenced above, teaches a device for separating a mixture of polynucleotides by ion pair-reversed phase-high performance liquid chromatography. The device comprises a polymeric monolith having non-polar chromatographic surfaces. The monolith comprises an underivatized poly-(styrene/divinylbenzene) matrix and is contained within a tube having an inner diameter in the range of 1 to 1000 micrometers.

4. Laboratory protocol notebooks regarding experiments related to this invention were kept by me as a Ph.D. student under the direction of Christian Huber.

5. I worked in Dr. Huber's laboratory during 1998 and 1999.

6. According to laboratory protocol notebooks submitted herewith, the first synthesis of PS/DVB monolith using decanol and tetrahydrofuran as porogens was performed on August 6, 1996. We then succeeded in a first separation of proteins (lysosome from beta-lactoglobulin B) in a PS/DVB monolithic column on August 25, 1998. See, e.g., the chromatograph at the bottom right-hand corner of the fourth laboratory notebook page. The first successful separation of oligonucleotides on a PS/DVB

BEST AVAILABLE COPY

Attorney Docket No.: P-576 (TI-0022)
Inventors: Huber et al.
Serial No.: 09/770,410
Filing Date: January 25, 2001
Page 3

monolith synthesized with decanol/THF as porogens was February 10, 1999.

7. We were able to fully practice our invention described in the above-referenced patent application prior to the date of the publication of the Gusev paper. A copy of the relevant laboratory notebook pages hereby accompanies my declaration.

I further declare that all statements made herein of my own knowledge are true and that all statements made on information and belief are believed to be true, and further, that these statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001 of Title 18 of the United States Code, and that such willful false statements may jeopardize the validity of the application or any patent issuing thereon.

24.02.2006

Date



Andreas Premstaller, Ph.D.

NOT AVAILABLE

6.8.98

Chloridkopolymer mit THF

C₁₀O₄H₁₆
mischkopolymer

Nr.	Datum	Kapillare ID/OD (µm)	Polymerisationsmischung					Temperatur [°C]
			Styrol (ml)	DVB (ml)	AIBN (g)	C12OH (ml)	THF (ml)	
M14_1	05.08.98	320/450 WBH06A, 20 cm, vs	1.00	1.00	0.050	3.00	0.00	70, TS
M14_2	05.08.98	320/450 WBH06A, 20 cm, vs	1.00	1.00	0.050	2.90	0.10	70, TS
M14_3	05.08.98	320/450 WBH06A, 20 cm, vs	1.00	1.00	0.050	2.80	0.20	70, TS
M14_4	05.08.98	320/450 WBH06A, 20 cm, vs	1.00	1.00	0.050	2.70	0.30	70, TS
M14_5	05.08.98	320/450 WBH06A, 20 cm, vs	1.00	1.00	0.050	2.60	0.40	70, TS

THF sollte kommen, in der Kette am Ende stehen als Toluol
THF distilliert, da mit Rückfluss (Theor.) stabilisiert

Ausgangsmaterial: VS 38.98 Tocken
THF dest.

Thet: 6.8.98 100h
Lk: 7.8.98 120h

T = 70°C

T = 70°C

BEST AVAILABLE COPY

7.8.98

M11_1 16 cm		
Fluß (µl/min)	Gegendruck (bar)	(bar/cm)
5	1	0.06
10	1	0.08
25	4	0.25
50	7	0.44
100	14	0.58
150	21	1.31
200	28	1.75
k (bar cm ³ µl ⁻¹ min ⁻¹)		0.008712

M11_2 15 cm		
Fluß (µl/min)	Gegendruck (bar)	(bar/cm)
10	1	0.07
25	4	0.27
50	8	0.53
100	14	0.93
150	20	1.33
200	25	1.67
k (bar cm ³ µl ⁻¹ min ⁻¹)		0.008303

M11_3 15 cm		
Fluß (µl/min)	Gegendruck (bar)	(bar/cm)
10	1	0.07
25	3	0.20
50	6	0.40
100	11	0.73
150	14	0.93
200	19	1.27
k (bar cm ³ µl ⁻¹ min ⁻¹)		0.008114

M11_4 16 cm		
Fluß (µl/min)	Gegendruck (bar)	(bar/cm)
5	3	0.56
10	14	0.88
25	32	2.00
50	68	4.25
100	125	7.58
150	150	11.25
k (bar cm ³ µl ⁻¹ min ⁻¹)		0.074902

M11_5 16 cm 200bar Gegendruck		
Fluß (µl/min)	Gegendruck (bar)	(bar/cm)
1	21	1.31
2	33	2.06
3	49	3.06
4	63	3.94
5	77	4.81
7	92	5.75
10	170	10.63
12	184	11.50
k (bar cm ³ µl ⁻¹ min ⁻¹)		0.963149

Ansat THF		Stärkung
% Porogen	k (bar cm ³ µl ⁻¹ min ⁻¹)	
0.0%	0.00871	
3.3%	0.00830	
6.7%	0.00811	
10.0%	0.07460	
13.3%	0.96315	

Ultrahoch:

55 Stufen: 450 cm

1 Ter → 8 µm

M11_1 kleine Poren ~ 3 µm
sehr porös.




M11_2 große Poren, nicht gleichmäßig

M11_3 große Poren, kleine Kanäle (µm-groß)

M11_4 sehr dicht, kleine Poren, porös, da nicht
ganz anders

M11_5 keine Ultrahoch zu erkennen.

el



es möchte sein, dass MA-4 und MA-5 genau stehen.
 für MA-4 und MA-5.

BEST AVAILABLE COPY

25.08.98

M 11.5 min 6.8.98

1 min 240 bar / 5 μ l/min

Fig. AP80875 S170
GYN 5000

SYNAP. 130 μ l/min \rightarrow Split \rightarrow 4.6 μ l/min
2 min 1000 bar / 10 μ l
2 min 30 sec 4 μ l/min

Equilibrium:

(A) H₂O 0.1% TFA

(B) ACN 0.1% TFA

50% A 14.50 -

Elution - T-Hold 2 min 15 sec T-Hold

10 μ l 2 min 15 sec $\frac{10}{2.75}$ 4.44 μ l/min

Baseline:

Thiohematop 0.05% [H₂O 50% ACN, 0.1% TFA

p = 200 bar

100% H₂O, 0.1% TFA: Proteine zeigen kein Peak \rightarrow kleine Kapillare?

Thiohematop near co. 1.50 min

50% ACN, 0.1% TFA: Proteine gleichzeitig mit Thiohematop: keine Reaktion
RTBA

27.8.98


50% ACN, 0.1% TFA: Protein quickly eluted, eluted as Peak.

LAC A.

LYS keine Reaktion

Now 40% ACN

BEST AVAILABLE COPY




50% ACN, 0.1% TFA
 of 2 ng/Lac Protein
 Retention on LacA bei 50% ACN

BEST AVAILABLE COPY

GP-Integration SYSL - C:\A980825.SMP Page 2
 950VB 100x0.32mm, C120H/TNF-Perogen, M11 3 98084 1998-08-27/20:13
 LysLacB 1mg/ml UV-1-1 1998-08-27
 Modified: 820/33-6.15 minACN/0.15TFA, 4.5/130ul/min, 25°C SynkeSoft V5.50
 Exp. No/Pos: 35/1 Control: Standard: 20.0 ul
 Sample Type: Integration Signals: AND11.SIG Inject: 1.00000
 Acquisition: 1998-08-27/19:53 Report: Weight: 1.00000
 Method: DEFAULT.LMT P-Table:


Effluente Proteinierung:
 Lys, LacB of 1 mg/L, 20ul 1/2.

30-60% ACN / 15 min, 0.1% TFA
 4.5 / 130 ul/min
 215 nm



No.	Ret. Time	Type	Area	Height	Half Width	Base Width	Plates
	min		mAU*min	mAU	min	min	
1	1.414	MS	3.351e+1	0.84	0.415	0.720	64
2	9.060	MS	3.170e+0	4.46	0.627	0.972	1156
3	13.379	MS	2.320e+0	4.82	0.567	0.988	1080
---	---	---	5.443e+0	10.12	---	---	---

A980825-36



Separation of proteins in a monolithic capillary column

Column, PS-DVB (monolith, 100 x 0.32 mm); chromatographic conditions, mobile phase, (A) H₂O, 0.1% TFA, (B) ACN, 0.1% TFA, linear gradient, 30-60% B in 15 min; flow rate, 4.5 µl min⁻¹; temperature, 25 °C; detection, UV, 220 nm; sample, lysozyme, β-lactoglobuline B, 20 ng each.

BEST AVAILABLE COPY

Hand line Trennung von Oligonucleotiden in Abwischen 1713.5

1713-5

$\lambda = 87 \text{ nm}$, $\text{id} = 200 \mu\text{m}$

Eluent: A: 50mM TEAA pH 6.8

B: 50mM TEAA 20% ACN pH 6.8

Temperatur: 50°C

Spaltkopf: TSP075375, 6cm Fluss 120/3.3/11min / 946

AG: A990209...ST11

Trennung von dT_8 , dT_{16}
 Testzeit mit Gradient 0-100% B/10min. 0.11 min
 = 6.65

Trennung mit dT_{12-18}


verschiedene Gradienten versucht.

gute Trennung: 30-50% B/10min

6-10% ACN/10min

Operator: c72551 Timebase: A990209 Sequence: A990209 Page 4-1
 10.2.1999 2:35 PM


4 dT80.24dT180.48ug/dl
 0-100%B/10min; A: 50mM TEAA pH 6.8; B: 50mM TEAA 20% ACN pH 6.8; 120/3.3/11min; D: 2min; 50°C
 Sample Name: dT80.24dT180.48ug/dl Injection Volume: 20.0 μL
 Control Program: Channel: UV-1-1
 Quantif. Method: OLIGO1 Recording Time: 09.02.99 19:00



No.	Ret. Time min	Area mAU*min	Height mAU	Half Width min	Plates (EP)	Asymmetry (AIA)
1	1.5242	0.450	4.053	0.165	13593	1.303
2	2.8154	0.744	8.025	0.111	10926	1.331
Total:		1.194	10.081			

Operator: c72551 Timebase: A990209 Sequence: A990209 Page 11-1
 10.2.1999 2:34 PM

11 dT12-18 0.25ug/dl
 30-50%B/10min; A: 50mM TEAA pH 6.8; B: 50mM TEAA 20% ACN pH 6.8; 120/3.3/11min; D: 2min; 50°C
 Sample Name: dT12-18 0.25ug/dl Injection Volume: 20.0 μL
 Control Program: Channel: UV-1-1
 Quantif. Method: OLIGO1 Recording Time: 09.02.99 21:29



No.	Ret. Time min	Area mAU*min	Height mAU	Half Width min	Plates (EP)	Asymmetry (AIA)
1	4.156	0.075	0.405	0.175	3142	1.050
2	4.707	0.071	0.384	0.178	3860	1.551
3	5.224	0.088	0.420	0.182	4101	1.248
4	5.709	0.073	0.370	0.180	5552	1.084
5	6.122	0.082	0.412	0.189	5813	n.a.
6	6.483	0.086	0.441	0.182	7042	n.a.
7	6.835	0.082	0.454	0.171	8586	n.a.
Total:		0.556	2.988			